

Autoclave Report

CE Softwall Cleaner™

International Article Number

162232

Product classification

Above the Floor Tool

Composition

Bushings:	Acetal
Composite Bolt and Nut:	Ultem
Inner and Outer Handle:	Aluminium
Frames & Hardware:	Stainless Steel



Test Conditions

Autoclave

- MLS Series 2.6 cu.ft. Portable Top-Loading Autoclave MLS-3781L-PA

Test Parameters

- 134 °C for 5 min
- Pre-vacuum: 20 kPa
- Cooling-down of product between autoclaving cycles

RESULT

After 10 autoclave cycles the tool does not show any damages or loss of function, but does show some characteristic changes: (See Below)

Label:

The label gradually faded over the autoclave cycles.



0 Cycles



10 Cycles

Bushings:

The bushings settled at 3% shrinkage after 7 cycles. The performance of the tool was not affected.

Date: 11/30/2018

All information contained herein is given according to our specifications and best knowledge without any warranty or guarantee and can be changed without further notice. The suitability of our products to our customers specific applications and conditions of use has to be determined by our customers. In particular, product users shall not be released from their duty to check all health, safety and environment relevant properties of the delivered goods under their specific conditions of use.